

STATIC PROGRESSIVE WRIST ORTHOSIS



FLEXION/EXTENSION

The wrist orthosis utilizes the “Mackie Hinge” for use in the management of soft tissue flexion or extension contractures. The joint has a removable knob for a lower profile fit and better patient compliance. The wrist orthosis is available in two pre-fabricated sizes (small and medium), custom to measurements, and custom to cast.

PRONATION/SUPINATION

The wrist orthosis utilizes the “Mackie Hinge” for use in the management of soft tissue pronation and supination contractures of the hand. The joint has a removable knob for a lower profile fit and better patient compliance. The orthosis is adjustable to fit a wide variety of patient sizes and is available in one prefabricated model only.

- Bi-directional infinite angle settings in fractional increments
- Reduce chances of micro-tears to soft tissue and the inflammation they cause
- Increase pronation and supination range of motion
- Patient friendly for easy home use
- Lightweight with removable knobs

**CALL US FOR ORDERING OR PRICING INFORMATION.
(877) 536.6107**